


Nutrient Pad Sets

Dehydrated Media Pads In Petri Dishes with Matching Membrane Filters



제약, 물, 화장품, 식품, 음료 산업 등에서 미생물 배양을 위해 사용하는 Petri Dish 형태의 배지와 격자무늬의 Membrane Filter가 함께 판매되는 제품

- 편리성: 미생물 Colony Counting을 위해 사용 하는 배지와 Membrane Filter가 함께 제공
- 멸균 및 건조 된 형태의 Nutrient Media로 장기간 (18개월 이상) 보관 가능
- Nutrient Pad위에 Membrane Filter를 올려놓고 배양하거나 샘플량이 적은 경우, 직접 도포 후 배양 가능




Ordering Information

Nutrient Pad Sets for Total Colony Counting

Order No.	Determination of	NPS Type (Filter Type)*
14063-47-N	Total count	Caso ⁽¹⁾
14084-47-RDN	Total count	R2A ⁽¹⁾
14084-47-N	Total count	R2A ⁽¹⁾
14055-47-RDN	Total count	Standard TTC ⁽¹⁾
14055-47-N	Total count	Standard TTC ⁽¹⁾
14085-47-N	Total count	Standard TTC I mod. ⁽¹⁾
14064-47-N	Total count	Standard ⁽¹⁾
14076-47-RDN	Total count	TGE ⁽¹⁾ Tryptone Glucose Extract
14076-47-N	Total count	TGE ⁽¹⁾ Tryptone Glucose Extract
14090-47-N	Total count	Yeast Extract ⁽¹⁾

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters (order no. -RDN = Microsart® e.motion Membranes)



Nutrient Pad Sets for E. coli, Coliforms and Enterobacteria

Order No.	Determination of	NPS Type (Filter Type)*
14087-47-RDN	E. coli and coliforms	Chromocult ⁽⁷⁾
14087-47-N	E. coli and coliforms	Chromocult ⁽⁷⁾
14082-47-N	E. coli	ECD ⁽²⁾
14053-47-RDN	E. coli and coliforms	Endo ⁽⁹⁾
14053-47-N	E. coli and coliforms	Endo ⁽⁹⁾
14097-47-N	Enterobacteria, E. coli	MacConkey ⁽²⁾
14068-47-N	E. coli and coliforms	m FC ⁽²⁾
14068-50-PDN	E. coli and coliforms	m FC in closed petri dishes ⁽²⁾
14067-47-RDN	E. coli and coliforms	Teepol Lauryl Sulphate ⁽²⁾
14067-47-N	E. coli and coliforms	Teepol Lauryl Sulphate ⁽²⁾
14056-47-RDN	E. coli and coliforms	Tergitol TTC ⁽²⁾
14056-47-N	E. coli and coliforms	Tergitol TTC ⁽²⁾

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters (order no. -RDN = Microsart® e.motion Membranes)

Nutrient Pad Sets Starter Kit

Order No.	Determination of	NPS Type (Filter Type)*
14095-47-N	E. coli and coliforms, total count, yeasts and molds	Mixed types: Endo, Standard, Wort ^(1, 2, 3)

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters

Nutrient Pad Sets for Other Faecal Bacteria

Order No.	Determination of NPS	Type (Filter Type)*
14051-47-RDN	Enterococci	Azide KF Strep(1)
14051-47-N	Enterococci	Azide KF Strep(1)
14057-47-N	Salmonellae	Bismuth Sulfite(1)

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters (order no. -RDN = Microsart® e.motion Membranes)

Nutrient Pad Sets for Non-Faecal, Pathogenic Bacteria

Order No.	Determination of NPS	Type (Filter Type)*
14075-47-RDN	Pseudomonas	aeruginosa Cetrimide(2)
14075-47-N	Pseudomonas	aeruginosa Cetrimide(2)
14074-47-N	Staphylococci, Staph.	aureus Chapman(2)

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters (order no. -RDN = Microsart® e.motion Membranes)

Nutrient Pad Sets for Yeasts and Molds

Order No.	Determination of NPS	Type (Filter Type)*
14061-47-N	Wild yeasts	Lysine(3)
14086-47-CCN	Yeasts and molds	Malt Extract(8)
14086-47-N	Yeasts and molds	Malt Extract(6)
14069-47-N	Yeasts and molds	Sabouraud(10)
14070-47-N	Yeasts and molds	Schaufus Pottinger m green yeast and mold(4)
14072-47-N	Yeasts and molds	Schaufus Pottinger m green yeast and mold(5)
14080-47-RDN	Yeasts and molds	Schaufus Pottinger m green yeast and mold(6)
14080-47-N	Yeasts and molds	Schaufus Pottinger m green yeast and mold(6)
14083-47-N	Yeasts and molds	Schaufus Pottinger m green yeast and mold(3)
14091-47-RDN	Yeasts and molds	Schaufus Pottinger m green yeast and mold(8)
14091-47-N	Yeasts and molds	Schaufus Pottinger m green yeast and mold(8)
14089-47-N	Yeasts and molds and bacteria	Wallerstein Nutrient WL Nutrient(2)
14058-47-RDN	Yeasts and molds	Wort(3)
14058-47-N	Yeasts and molds	Wort(3)
14092-47-RDN	Yeasts and molds	Wort(8)

individually, sterile packaged in petri dishes, 100 per box, with 100 individually, sterile packaged 47 mm membrane filters (order no. -RDN = Microsart® e.motion Membranes)

* Filter Type

- (1) Green with dark-green grid, 0.45 µm pore size
- (2) White with green grid, 0.45 µm pore size
- (3) Gray (after wetting black) with white grid, 0.65 µm pore size
- (4) White with green grid, 0.65 µm pore size
- (5) White with green grid, 1.2 µm pore size
- (6) Gray (after wetting black) with white grid, 0.8 µm pore size
- (7) White with black grid, 0.45 µm pore size
- (8) Gray (after wetting black) with white grid, 0.45 µm pore size
- (9) White with green grid, 0.45 µm pore size, High Flow
- (10) Gray (after wetting black) with white grid, 0.45 µm pore size, High Flow